**A new lens for comfortable near and intermediate vision.**

**Rodenstock Nexyma**

### Content

1. Product details
2. Product comparison
3. General hints for ordering
4. Lens calculation / optimization
5. Production
6. Complaints
7. Enlargeable delivery range

- Nexyma offers the chance to organize our near lenses in a new way.

**Nexyma 80:** lens for intermediate and near vision (successor of Office)

**Nexyma 40:** reading lens (successor of Cosmolit P)

- One product name
  - new category „near vision lenses“
  - One simple ordering method
  - All near lenses have a clear near vision.
- Transparent and easy to use system.

# Nexyma Marketing Concept

## Starting position

---

### Product description

Introduction of two different near-vision lenses which cover different main vision areas:

- Nexyma 40  
for a main vision distance of 40 cm up to approx. computer monitor distance
- Nexyma 80  
for a main vision distance of 80 cm up to near room distance


# Nexyma Marketing Concept.

## Communications strategy

---

### Target group

Presbyopes from around the age of 40 who work a lot on activities in the near and intermediate area.

# **Nexyma Marketing Concept.**

## **Communications strategy**

---

### **Core message.**

Everything in the extended near area can be seen in optimum width and clarity with a relaxed head and body posture with Nexyma near-vision lenses.

### **Reason why.**

The vision zones of Nexyma near-vision lenses have been optimised for the natural lowering of the gaze for all activities in the near/intermediate area (40cm to approx. 120 cm). Reading glasses and normal progressive lenses are subject to restrictions here.

Advantages for the optician	Advantages for the consumer
<ul style="list-style-type: none"><li>• easy ordering</li><li>• transparent system for near vision lenses</li><li>• easy lens selection by application</li><li>• separate price list category: easy to find</li><li>• good chance to sell a 2nd pair of lenses</li></ul>	<ul style="list-style-type: none"><li>• Rodenstock offers the right lens for individual vision habits in the near and intermediate distance</li><li>• comfortable half eye frames possible</li><li>• higher vision comfort than with normal SV reading lenses offer</li><li>• fatigue free vision</li><li>• relaxed head posture</li></ul>

**The right lens for all visual demands for near and intermediate distances dependent on personal preferences.**

Communication to optician	Communication to consumer (printing data available)
<ul style="list-style-type: none"><li>• trade advertising</li><li>• product folder</li><li>• centring card</li><li>• advice card</li><li>• additional card for demo viewer</li><li>• separate price list category</li></ul>	<ul style="list-style-type: none"><li>• direct mailing</li><li>• consumer leaflet</li><li>• POS material</li><li>• poster</li></ul>

# Nexyma

## Material and range

---

	Nexyma 40	Nexyma 80
Material	Perfalit 1.5	
Refractive index	1,502	
Abbe No	58,2	
Density	1,32 g/cm <sup>3</sup>	

### Power Range (strongest meridian)

Ø 65	-4.00 D up to +8.00 D
Ø 70	-4.00 D up to +6.50 D
Ø 75	-4.00 D up to +5.00 D
Cylinder	up to +4,00 D
Prism	up to 5 cm/m

# Nexyma

## Technical Data

	Nexyma 40	Nexyma 80
Optimized for main vision distance	40 cm	80cm
Ordering	With near power at 40cm	
Degression up to add 1.75 D	1.0 D Engraving N4	0.8 D Engraving N8A
Degression from add 2.00 D		1,50 D Engraving N8B
Min. fitting height	13	20
Min. frame height	19	26
Optimized for the position of wear	yes	
Tints and coating	As available for Perfalit 1.5	

### **The fitting and ordering is as easy as always:**

Centre the centring cross to the middle of the pupil centre with a habitual head and body posture and the eyes looking straight ahead to the distance.

For the order just specify the near refraction at 40cm and the lens type:

Nexyma 40	from add 1.25 D	age from approx. 45 years
-----------	-----------------	---------------------------

Nexyma 80 A	up to add. 1.75 D	age from 45 to approx. 50 years
-------------	-------------------	---------------------------------

Nexyma 80 B	from add. 2,00 D	age from approx. 50 years
-------------	------------------	---------------------------


# Nexyma Marketing Concept.

## Communications strategy

---


### Progressive lens

Natural head posture for distance vision mainly  
Less comfortable for extended tasks in the near and intermediate vision area


### Nexyma 80

Optimum solution for the vision in approx. 80 cm  
Fatigue free working due to natural head posture  
Additional clear near vision


### Nexyma 40

Optimum solution for the vision in approx. 40 cm  
Ideal for relaxed long term reading  
Significant enlarged vision area over the near area only

## Stamp Markings/Engravings semi finished lenses:


**Engraving: horizontal ellipse**

„N4“ for Nexyma 40

„N8“ for Nexyma 80

## 1. Stamp marking and engraving finished lenses

- The function engraving is a flat ellipsoid.
- The engravings are positioned symmetrical to the centring point and at the same height as the horizontal line of the lens.
- PRP is placed 1mm nasally, (lens is not decentred)
- Fitting like progressive lenses
- NRP for **Nexyma 40**: 7mm below the horizontal line of the lens
- NRP for **Nexyma 80**: 14mm below the horizontal line of the lens


### 1. Nexyma 40

- Min. fitting height: 13mm
- Min. frame height 19mm for full use of the degression
- Nexyma 40 can be used with half eye frames  
For more narrow frames (<19mm) the centring cross should remain within the lens. Usually that results in degression loss.


### 2. Nexyma 80

- Min. fitting height: 20mm
- Min. frame height 26mm for full use of the degression.  
Nexyma 80 is recommended for full frames only.


**Patient: age < 50 years (add 1.50)**

**Conv. PAL**

**Nexyma 80A**


**Nexyma 40**


**Patient: age > 50 years (add 2.50)**


**Conv. PAL**

∞ ☺☺☺

**Nexyma 80B**

∞ ☹


**Nexyma 40**


∞ ☹


Conventional progressive lens  
Add 1.25

Nexyma 40


Nexyma 80 A

The lower the blue intensity is the lower the refractive errors and the larger the useful vision areas are for a certain distance.

## 1. Product details

Sabine:  
Multigressiv 2  
sph +0.50 Add 1.25


Black numbers:  
(comfort area)

Blue numbers:

Distance using „normal“ accommodation


Distance using maximum accommodation

Franz:  
Multigressiv 2  
sph +0.50 Add 2.50


## 1. Product details

Sabine:  
Nexyma 40  
Nahwert +1.75D

Franz:  
Nexyma 40  
Nahwert +3.00D

Black numbers:  
(comfort area)

Blue numbers:


Distance using „normal“ accommodation

Distance using maximum accommodation

## 1. Product details

The degression for **Nexyma 40** is 1,00D.


The necessary near refraction will always be reached. The size of the vision area results as shown in below graph dependent on the rest accommodation (addition).

 Main vision area  Enlarged vision area

## 1. Product details

Sabine:  
Nexyma 80 A  
Nahwert +1.75D


Black numbers:  
(comfort area)

Blue numbers:

Distance using „normal“ accommodation

Distance using maximum accommodation

Franz:  
Nexyma 80 B  
Nahwert +3.00D


In order to reach the near refraction power in 40cm and at the same time also the refraction power for the main vision area in 80cm it is necessary to have two degression versions.

Dependent on the addition an enlarged vision area is reached.


Cosmolit

P


Due to the high peripheral refractive errors Cosmolit P is not up to date anymore.

Nexyma 40

⇒ lower peripheral refractive errors  
⇒ larger depth of vision  
⇒ significant wider clear near vision areas

Nexyma 80 B

⇒ universal usable  
⇒ 100% near power


# Nexyma 40

technical product information

## Power gradient Nexyma 40 Degression 1.0 D

### Reference points

- 1 Reference point near  $B_N$ =measuring point for optician
- 2 Lens horizontal
- 3 Centering Cross  $B_Z$


# Nexyma 80A

## technical product information

### Reference points

- 1 Reference point near  $B_N$ =measuring point for optician
- 2 Lens horizontal
- 3 Centering Cross  $B_Z$

### Power gradient Nexyma 80 A Degression 0,8 D up to Add 1,75D


# Nexyma 80B

## technical product information

### Reference points

- 1 Reference point near  $B_N$ =measuring point for optician
- 2 Lens horizontal
- 3 Centering Cross  $B_Z$


### Power gradient Nexyma 80 B Degression 1,5 D from Add 2,00 D


## 2. Product comparison

How can you differentiate different products?

- a) in design
- b) in power increase respectively length of increase (digression)
- c) resulting lowering of gaze for the near vision

## 2. Product comparison

Product	Design	Power increase [D] at a length [mm]	Near power reached at:
Zeiss Business	C A	1 D / 30mm (30mm)	Approx. y= - 14mm
Zeiss Gradal RD	C A	2 D / 30mm (30mm)	Approx. y= - 17mm
Essilor Interview	C D C D	0.8 D / 20mm (12mm) 1.4 D / 20mm (12mm)	Approx. y= - 10mm Approx. y= - 10mm
Nexyma 40	D	1.0 D / 15mm (15mm)	y= - 7mm
Nexyma 80 A 80 B	B	0.8 D / 22mm (22mm) 1.3 D / 22mm (22mm)	y= - 14mm y= - 14mm

## **2. Product comparison**

Difference compared to conventional progressive lenses:

- lower peripheral refractive errors
- larger vision zones in the main vision area.

For each Nexyma version this results in a significant advantage over progressive lenses in the certain field of application .

Difference to other near vision lenses:

- Advantages in the width of the visual fields.
- Nexyma offers an adjusted lowering of the gaze for the main vision task
- Comfortable and non tiring usability of the lens possible.

### **3. General hints to the order procedure**

- Orders are placed with the near power,
- Data entry like “normal distance vision” Order with distance + addition
- If there is an addition specified - Nexyma is ordered like a progressive lens, this addition value is ignored when the lens is calculated.
- EDV Code:  
NEV = Nexyma 40  
NEA = Nexyma 80 A  
NEB = Nexyma 80 B
- SF lenses: ordered by product code and BC
- Single lens orders are possible like with progressive lenses.

## **5. Peculiarity for the production**

No peculiarity.

Can be produced like e.g. Progressiv life2.

Limitations can occur dependent on the place of production.  
(Keyword: Blocking rings etc.)

## **6. Complaints / Rejects**

Intolerance because of missing near power can be eliminated when the refraction data was determined correctly.

When the vision areas are insufficient or the depth of vision is too low, please check:

- Does the working environment correspond to the profile of the product 40 or 80?
- Is the patient able to use the main vision areas in the right height?
- Did you choose the right type of degression for Nexyma 80?
- Is the power for near vision correct (too high/not necessary)?

## **7. Extendable power range**

- No near powers higher than +8D and smaller than -4D possible (With BC pre-setting extension might be possible, Limitations have to be expected due to the “BC effect” especially with conv. progressive lenses – therefore a previous judgment is necessary).
- No higher cylinders than 4 D possible
- In cases of anisometropia the base curves R/L are not adjusted.

**Focus: The benefit for the consumer.**

### **1. Realised using the headline:**

**Take a more relaxed view of life!**

**Nexyma near-vision lenses go further than reading glasses.**

- Optimum clarity of vision in the near and intermediate area.
- With a relaxed head and body posture
- Ideal for wear at home, for hobbies and at work.

# Nexyma Marketing Concept.

## Communications strategy


---

**Focus: The benefit for the consumer.**

### **2. Realised using new image material:**

Shows the deep breadth of different possibilities for use

Visualisation of the benefits in different situations  
(E.g.: reading, working at a desk, hobby)

# Nexyma Marketing Concept.

## Communications strategy

---

**Focus: The benefit for the consumer.**

### **3. Realisation using new demonstration tools:**

New slide card for the 3D Demo-Viewer


# Scene with natural vision

# Vision with progressive lenses

# Vision with Nexyma 80

# Vision with Nexyma 40

**R**  
RODENSTOCK

# Nexyma Marketing Concept.

## Communication tools for the optician

---

### Mailing of the folder

Contents:

**Optician's  
Sales folder DIN A4**

# Nexyma Marketing Concept.

## Communication tools for the optician

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